

## Manual for the "Biologist" category

### Table of contents

Foreword .....	2
Home page .....	2
Emergency management .....	3
Biologist work list .....	3
Group validation .....	4
Validation by file .....	6
Status of current withdrawals .....	6
List of files .....	7
System settings .....	8
Benchmark of analyses .....	9
Use of the repository via the software .....	9
Add a new analysis .....	10
Import of repository .....	10
User management .....	11
Exporting users and connections .....	11
Reports .....	12
Epidemiological Report .....	12
Statistical report .....	13
Patient history .....	14
End of day report .....	15
Daily billing statement .....	15
Quality and Non-Compliance Module .....	16
APPENDIX: Corresponding ID labels - Data export .....	17

## Foreword


This manual presents the elements of LabBook accessible to a person with "biologist" rights. If you do not have access to any of the actions via your interface, please contact your administrator so that these rights can be assigned to you.

For a biologist's work to be finalized for LabBook, the analysis result must be validated biologically. Once this validation is done, a report appears to provide feedback to the patient.

## Home page


The screenshot shows the LabBook home page interface. At the top, there is a navigation bar with icons for Administrative, Reports, Settings, Quality, and Non-conformity. The user's name 'Bernard BIO' and a language selector 'US' are visible in the top right. Below the navigation bar, there are several main sections: 'Record list', 'Biologist work list', 'Status of specimens in progress', and 'Results entry'. A callout box points to the user profile area, stating 'Allows you to disconnect' and 'Edit profile and password'. Another callout points to the language selector, stating 'Allows you to change the language of the'. A third callout points to a backup icon, stating 'Here you have the information about the last backup'. Below the navigation bar, there is a pie chart titled 'Validations to be carried out' and a section titled 'Laboratory' with a list of items.

By hovering the mouse cursor over the icon next to your name, you have the option to log out, edit your profile and password:




- Disconnection: Allows you to leave the application by closing your account.
- Edit user: Allows you to edit your profile information (login, first name, last name, ...)
- Edit password: Here you can change your password

This home page allows you to directly access the file entry and work lists for your profile.





The "Overall Report" submenu allows you to download all the reports in a single pdf file. You have the possibility to specify a date and also "Exclude those already downloaded".



For the use of the functions "Results entry" and "Current sample status", please refer to the Manual for the "Technician" category.

## Emergency management

If you see the logo  on your home screen, it means that some files have analyses to be done urgently. By moving the mouse over this logo, the software will indicate the number of analyses



concerned:

You have the list of these analyses in urgency by clicking on the icon "Urgencies".

## Biologist work list

(In the "Administrative" menu or directly on the "Biologist work list" icon on the home page)

This list shows you the files that require your intervention. This list allows you to make a grouped biological validation by clicking on the button "Start the grouped validation" : that is to say that the

files will be linked to each other to allow you to validate them directly one after the other.

**Work list**

Search

Record date greater than   Record date less than   Urgent only

Type  Patient code

**Record list**

Total number of lines : 2  
First Previous 1 Next Last

Action	Status	Emergency	Record number	Record creation	Patient code	Name	Firstname
			0007	2023-08-08	TOHU4	Fall	Tapha
			0006	2023-08-07	ELIB5	KA	OUMAR

**Group validation**

## Group validation















































































The grouped validation allows you to link the validations for the files requiring an action from you. You will therefore enter each file to validate it. Once this is done, you will move on to the next file. At the end of all the validations, a grouped printout will be possible using the button :

(Which can be found at the bottom left of the page)

Validation biologique - 169 QAWG6 - test RAKOTO - Né(e) le 2000-03-21 - 22 années

Hémogramme (Numération Formule Sanguine) [Hématologie, Immunohématologie et Hémostase]

03/2023

Numération			
			
			
			
			
			
			
			
			
			
			
Formule leucocytaire			
			
			
			
			
			
			
			
			
			
			
			
			
			
			
			
			

Once the validation is done, you can :

- Add new analyses if you think additional analyses are needed, by the option at the bottom of the page;


Add new analyzes

Click to start a search

Next record


- Go to the next folder with the button (at the end of the page and on the right)

The report is generated on each page at this location:



And by clicking on this link, you have the possibility to download and print the report. You can see the number of downloads of the report with the number (..) after the Model result.

## Validation by file



For each folder, in the "Action" column, click on the down arrow next to the "i" icon, you will have several choices:

- "Administrative file": allows you to view the information entered when the file was registered. It also allows you to print an invoice or to reprint a report.
- "Biological validation": allows the biological validation of the recorded results and the output of an associated report.
- "Delete" : BE CAREFUL ! This allows you to delete a folder.
- The biologist account has the possibility to do the tasks of the secretary account through the "Results entry" action and the tasks of the technician through the "Technical validation" action.

## Status of current withdrawals

This module allows you to enter the samples taken in the laboratory. To do so, you can go through the icon on your home page or through "Administrative", "Status of specimens in progress".

Record number	Record date	Patient	Blood	Stool	Urine	Other
0007	2023-08-08	Fall Tapha	To complete			
0006	2023-08-07	KA OUMAR	To complete			
0005	2023-08-05	ONE Trois	To complete			
0004	2023-08-05	TWO Deux				To complete

This table summarizes for each file, the list of withdrawals to be made.

If the line shows "To be entered", it means that no data has been entered for the collection in question.

If the line shows "Modify", it means that information has been entered for the file. However, you can access it at any time to post a change.

## List of files

You can access this list either via the icon on the home page or via the menu "Administrative", "Record list".

**Record list**

Name  Firstname  Code  Status

Type  Urgent only

**Record list**

Total number of lines : 16  
First Previous 1 Next Last

Action	Status	Emergency	Record number	Record creation	Patient code	Lab code	Name	Firstname
	B		0007	2023-08-08	TOHU4	12	Fall	Tapha
	T		0006	2023-08-07	ELIB5	14211-cares	KA	OUMAR
	B		0005	2023-08-05	KFQY5		ONE	Trois
	B		0004	2023-08-05	ZIYT7		TWO	Deux

This display allows you to quickly see the files in progress. So if someone calls you to know the status of a file you can quickly answer. The different existing statuses are the following:


- A: administratively validated (no test results have been entered in the file)
- T: technically validated (results have been entered and validated by the technician. The biologist must validate the results to create the report).
- B: biologically validated (the results were validated. The report could be edited. By clicking on the file, you can see when).
- I: some elements have been validated but not all of the requested analyses.

N.B.: If your file has an "Internal laboratory file number", it is displayed in the file number column

(and after the file number)

In the Actions menu, on , you will find the elements that allow you to act on the file as presented in the technician work list.

## System settings



Via the menu "Settings", "Preferences", you can make changes to the parameters of your application. Each line can be edited directly on the input fields of the value column.

Everything about document headers allows you to change the way your lab is displayed on invoices and reports.

The screenshot shows the 'List of preferences' settings page. The page has a purple header with the title 'List of preferences'. Below the header is a table with two columns: 'Label' and 'Value'. The table contains various settings, including document headers, inpatient billing, and reporting options. At the bottom of the page, there are 'Exit' and 'Save' buttons.

Label	Value
Unit price for sampling and analysis procedures	1000
Document header 1	Nom du laboratoire
Document header 2	Sxxx au capital de xxx € RCS xxx xxx xxx autorisation n°xx xxx
Document header 3	Horaires : du lundi au vendredi : 07h00-19h00, le samedi : 07h30-12h00 www.example.com
Document header - Address	100 place de la République 10000 Maville
Document header - Phone	01 23 45 67 89
Document header - Fax	01 98 76 54 32
Document header - Email	labo@examples.com
Document header - City	Ma ville
Inpatient billing	<input checked="" type="checkbox"/>
Default age unit (Days, Weeks, Months, Years)	Years
Automatic disconnection (mn)	23
Quality module activated	<input checked="" type="checkbox"/>
Invoicing module activated	<input checked="" type="checkbox"/>
Language of reports	French
Language of the repository	French



If your laboratory's practice is to bill for tests performed on hospitalized patients, you only need to edit the line corresponding to "Inpatient billing" to replace the 0 with a 1.

The Quality and Invoicing modules can be activated/deactivated in the same way by replacing the 0 (inactive) with 1 (active).

The automatic logout time allows you to determine after how many minutes of inactivity a user is logged out. This ensures that a user does not stay logged in for too long when not using the software, which could allow someone else to access the software on their behalf.

You can also change the languages of the reports and the repository. For the repositories, there are only two languages available: French and English.

The "Save" button allows you to save these settings.

## Benchmark of analyses

The test repository is the element that allows you to indicate the tests done in your laboratory. It allows you to activate/add/deactivate them. It also allows you to change their prices, their analysis family and the sample(s) needed for this analysis.

### Use of the repository via the software

Designation of the act

Analysis family

Type of specimen

Active analysis: Yes

Total number of lines : 508  
First Previous 1 Next Last

Action	Code	Designation	Abbreviation	Family	Status	Bio. product
	781	Dépistage syphilis		Immunologie	Activated	
	ABCL	Antibiogramme 1ère ligne des mycobactéries en milieu liquide	ATBBKML TUB	Bactériologie	Activated	
	B001	Acide urique (uricémie)		Biochimie sanguine	Activated	PB1 : Prélèvement de sang veineux
	B002	Bicarbonates		Biochimie	Activated	PB1 : Prélèvement de sang

Search

At the bottom of the list, you have the possibility to add an analysis. Remember to check via the search engine that the analysis is existing or not before adding it so as not to create confusion when creating the folder.

First Previous 1 Next Last

Back Import the repository Export the repository List of variables Add an analysis

N.B.: The "List of variables" button allows you to see all the existing variables in the software and identify the number of analyses using each variable

## Add a new analysis

The screenshot shows the 'Add a new analysis' form. At the top, there is a navigation bar with icons for Administrative, Reports, Settings, Quality, and Non-conformity. The user profile 'Bernard BIO' is visible in the top right. The form itself has a purple header 'Analysis' and a sub-header 'Analysis'. The fields include: Code \* (text input), Designation of the act \* (text input), Abbreviation (text input), Analysis family (dropdown), Type of specimen (dropdown), Rating unit (text input), Quotation value (text input), Active analysis (radio buttons for Yes and No), Whonet export (radio buttons for Yes and No), and a large text area for Comments.

The interface allows you to create new analyses by filling in at least the "Code" and "Procedure designation" fields.

You are also asked to indicate the variables that people should fill in for this analysis.

The screenshot shows the 'Variables' form. It features a search bar for variables with a dropdown and a green '+' button. The form includes fields for: Label \* (text input), Var. code (text input), Id (text input), Result type \* (dropdown), Description (text input), Normal value min. (text input), Normal value max. (text input), Underline (radio buttons for Yes and No), Formula (text input), Unit (dropdown), Accuracy (text input), Num. var for the formula (text input), Display position (text input), Comments (text area), and Mandatory result (radio buttons for Yes and No). There is also a 'Generate a QR code' option with radio buttons for Yes and No. Below the form is a table with columns: Action, Name, Unit, Min, Max, Num. var, and Position. The table is currently empty. At the bottom, there is a 'Cancel' button and a 'Save' button. The footer shows 'Version : 3.3.11 | Contributors'.

## Import of repository

First Previous 1 Next Last

Back

Import the repository

Export the repository

List of variables

Add an analysis

At the bottom of the repository list, you can import a repository. This allows you to load a repository provided on the LabBook site.

You can also export your repository (CSV format) by clicking on the "Export repository" button and modify the repository with a text editor and then re-import it.

## User management

Action	Id	Login	Firstname	Name	Status	Creation date	Role	Department	Origin
1	1	root	Abdel Assim Dembelé	root	Activated	2015-12-15	Administrateur		
2	2	biologiste	Bernard	BIO	Activated	2021-03-04	Biologiste		
3	3	technicien	Thierry	TECH	Activated	2021-03-04	Technicien		
4	4	techav	Thomas	TECHAVANCE	Activated	2021-03-04	Technicien avancé		
5	5	techq	Thibault	TECHQUALIT	Activated	2021-03-04	Technicien qualitticien		
6	6	secrtaire	Sophie	SECR	Activated	2021-03-04	Secrétaire		


This interface allows you to search for a user. In the actions menu, you can :

- Edit a user: this allows you to change the information (name, first name, email address) as well as its rights;
- Edit password: allows you to manually change the password of a user;
- Deactivate user: Allows you to suspend a person's access;

## Exporting users and connections


You can export the users to a csv file by clicking the "Export users" button at the bottom right of the list.

The "Export connections" button allows you to obtain a list, in a csv file, of the connections of each user. You can also specify the recovery of data between a date.










## Reports

You have access to the tool that presents the standard LabBook reports.




## Epidemiological Report

This report presents a number of indicators. You have the details of the report in the Epidemiology tab with the disease, the nature of the sample, the total received, the total analysis and the results.

Lab | Book  Administrative  Reports  Settings  Quality  Non-conformity Bernard BIO  Bio 

**Epidemiological report**








Date from  to   Search

## Epidemiology

Disease	Type of specimen	Total received	Total analysed	Results
Choléra	Selles	0	0	<b>Examen direct</b> Bacilles à mobilité polaire : 0 BGN incurvés : 0 <b>Culture</b> Négatif : 0 Positif : 0 V. cholerae O1 : 0 V. cholerae O139 : 0
Ebola	Sang	0	0	Positif : 0 Négatif : 0
Ebola	Crachat	0	0	Positif : 0


You can also search for a specific report by entering the date range on the "Search" tab.

## Statistical report

Lab | Book  Administrative  Reports  Settings  Quality  Non-conformity Bernard BIO  Bio 

**Statistic report**

Search

Date from  to  Requesting department   Search

## Distribution of records

	Less than 5 years	5 to 20 years	20 to 40 years	40 to 50 years	Greater than 50 years	No information	Total
<b>Male</b>	0	6	8	0	0	0	14
Female	1	0	1	0	0	0	2
<b>Unknown</b>	0	0	0	0	0	0	0
Internal	0	6	0	0	0	0	6
<b>External</b>	1	0	9	0	0	0	10
Total	1	6	9	0	0	0	16

In this report, you have the data concerning the distribution of the files by sex, age group, and type of patient.

The search for data is possible by date and also by requesting department for inpatients. The addition of the requesting services is done

You also have statistical figures on prescribers, samplers and samples. And at the end, a table that shows you the number of patients and analyses by type (outpatient, inpatient and on-call)

	M	F	U	Total		M	F	U	Total
Number of outpatients	13	2	0	15	Number of external scans	18	4	0	22
Number of hospitalized patients	1	0	0	1	Number of inpatient analyzes	3	0	0	3
Number of on-call patients	0	0	0	0	Number of Guard Scans	0	0	0	0
Number of patients	14	2	0	16	Number of analysis	21	4	0	25

## Patient history

This report allows you to search for a patient and access their history in the software.

**Patients list**

Search

Code  Name  Phone


Lab code  Firstname  Sex

Total number of lines : 9  
First Previous 1 Next Last

Action	Code	Lab code	Name	Firstname	Date of birth	Sex
<i>i</i> ▾	HFQF6					Unknown
<i>i</i> ▾	UGQL7					Unknown
<i>i</i> ▾	DHB11		AKAKPO		1990-07-03	Male
<i>i</i> ▾	TOHU4	12	Fall	Tapha		Male
<i>i</i> ▾	ELIB5	14211-cares	KA	OUMAR	1985-02-01	Male
<i>i</i> ▾	HMKA2		KOUDOUAZAN		2000-10-29	Female
<i>i</i> ▾	KJOO4		MANZI		2002-07-28	Male
<i>i</i> ▾	KFQY5		ONE	Trois		Male
<i>i</i> ▾	ZIYT7		TWO	Deux		Female

First Previous 1 Next Last

You can search for the patient concerned by his or her code (the code generated by LabBook), lab code (the internal patient code of the laboratory), first name and last name, telephone number and gender. To access the history, simply click on the down arrow next to the "i" icon in the Action



column . You will then arrive on a page presenting all the analyses done by the patient and the known results.

## Patient history

Code KFQY5  
 Lab code  
 Name ONE  
 Firstname Trois  
 Date of birth  
 Sex Male

### Analyzes list

Total number of lines : 12  
 First Previous 1 Next Last

Action	Record number	Prescr. date	Analysis	Variable	Result
	0005	2023-08-04	D Dimères par technique agglutination de particules de latex	D Dimères	12
	0003	2023-08-04	Acide urique (uricémie)	Acide urique	14
	0006	2023-07-26	Recherche de plasmodium (goutte épaisse et frottis mince)	Densité parasitaire	1.00

## End of day report

By default, this report shows the status of the various files in progress for the current day.

## End of day report

### Search

Date from 07/08/2023 to 08/08/2023 Search

Total number of lines : 3  
 First Previous 1 Next Last

Action	Date	Record number	Family	Analysis	Validation
	2023-08-07	0006	Immuno- Sérologie et Biologie moléculaire	Recherche des Ac anti-VIH (test rapide)	Administratif, Technique
	2023-08-08	0007	Biochimie sanguine	Glucose	Administratif, Technique, Biologique
	2023-08-08	0007	Biochimie sanguine	Urée	Administratif, Technique, Biologique

First Previous 1 Next Last

Exit

Export data

It is possible to change the filter on the dates to select another day or directly another period.

The "Export data" button allows you to obtain the report in csv format.

## Daily billing statement

This report allows you to view the billing status of cases.

Lab | Book Administrative Reports Settings Quality Non-conformity Bernard BIO Bio US

**Daily billing status**

Search

Date from 08/08/2023 to 08/08/2023 Record creator Click to start a search

List of current invoices Search

Total number of lines : 1  
First Previous 1 Next Last

Record number	Invoice no.	Price	To be paid after discounts
202308080001		20000.00	20000.00

First Previous 1 Next Last


Left to pay 0.00

Download the report

It is possible to modify the filter on dates as well as to restrict on a file creator in order to reduce the view by creator.

The outstanding balance is displayed at the bottom of the listing.

## Quality and Non-Compliance Module



By clicking on the icon the declaration of a non-conformity is possible for the authorized personnel.

The form that opens allows you to describe the type of non-conformity, its impact on the laboratory and the elements put in place to resolve this non-conformity.

For all other information on these two modules, please refer to the manual on these features.



## APPENDIX: Corresponding ID labels - Data export

<i>ID/Code</i>	<i>Wording</i>
<b>Sex (Gender)</b>	
1	Male
2	Female
3	Unknown
<b>impact_patient (Impact on the patient)</b>	
1053	Low
1055	Important
1057	Grave
0	No
<b>Impact_user (Impact on staff)</b>	
1053	Low
1055	Important
1057	Grave
0	No
<b>Title</b>	
260	Mr.
261	Ms.
262	Miss
263	Doctor
264	Professor
<b>Ctq_type_val (Control type)</b>	
QN	Quantitative
QL	Qualitative
<b>cte_conform (Control result)</b>	
C	Compliant

# Lab | Book

N	Not in compliance
U	Other